

Defeating Pain with MR4 Laser



Covers 5 times more area to treat large muscle groups
Incorporates 4 super pulsed GaAs diodes

Pulse drives the light energy, to the target tissue, up to 10-13cm (4-5 inches) deep

Revolutionary MR4 Super Pulsed Laser Technology for treatment of pain and injuries

Advanced Non-invasive Pain Relief Device

- Anti-inflammatory action
- Increases vascular activity/metabolic activity
- Stimulates nerve function
- Rapid cell growth
- Reduces fibrous tissue formation

Indication includes:

- Acute, sub-acute and chronic pain
- Shoulder pain
- Carpal tunnel
- Back pain, Tennis elbow
- Arthritic pain
- Many other indications

4 Radiances Technology

Combined into a **SINGLE** unit

- **The GaAs Super Pulse Laser / 905nm**
- **Four infrared diodes / 875nm**
- **Four red diodes / 660nm**
- **Static Magnetic Field/enhance energy potential at the molecular levels**

| Item # | Product Description | CAD Price |
|---------|--|---|
| 17-9211 | MR4 Laser (SE25) | \$4,295 → \$3,995 |
| 17-9212 | MR4 LS50 Shower Laser (4 x 905nm diodes) | \$5,995 → \$5,795 |
| 17-9213 | MR4 Laser Stim with Target Technology | \$6,795 → \$6,495 |
| 17-9214 | MR4 ULTRA Package Laser (Stim+Shower) | \$8,695 → \$8,395 |
| 17-9210 | TerraQuant TQ Solo | \$2,195 → \$1,995 |

Super Pulsed Laser Power

Radiation Penetrates
Peak Pulse Power 25,000mW
or 50,000mW in Laser Shower (20 cm²)
10-13cm (4-5 inches) deep

Also Great for Anti-Aging and Cosmetic Solutions

Advantages of MR4 and LaserStim™ with TARGET™

- Synergistic effect laser therapy provides benefits of light therapy and neurostimulation, identifying optimal treatment sites, and provides highly targeted treatment where photons need to be delivered and electro-stimulation needs to be applied.
- Faster pain relief and reduction in recovery time
- Constantly changing impulses, preventing accommodation of tissue to neurostimulation
- Deepest tissue penetration of photonic energy due to Super Pulsed Laser
- Preset frequencies allow clinician to control depth of therapeutic energy from the surface up to 13 cm (5 inches)
- Can also target acupuncture points

Technological Advantage: LaserSweep™

The continual use of the same frequency, as found in continuous wave lasers, can lead ultimately to the beginning of adaptation processes and thereby can diminish the therapeutic effect. The Multi Radiance MR4™ is capable of creating unique and exclusive, patented "sweep" laser programs.

With each laser impulse of a Sweep Program, the "frequency" is changed, varying between end points of pre-selected frequencies. This effectively change the laser's target depth of penetration, in essence, "scanning" the target tissue from superficial to deep, while preventing biological adaptation to the treatment.

ACCESSORIES FOR YOUR TERRAQUANT



17-9210 Portable TQ Solo

\$1,995

Laser Peak Power: 15W (15,000mw)

Infrared LED's: 60mW

Red LED's: 7.5mW

Static Magnetic Field

Wavelength: Laser 905nm,

Infrared 875nm, & Red 660nm

Class II Equipment Applied Part Type BF

Laser Class: 1M (IEC 60825-1)

Weight: 300g , Dimensions: 5.5 x 2 x 2"



Most versatile, perfect for use with probes

Provides complete spectral irradiance

17-9202 Extra SE25 Standard (TQ1) Emitter

\$1,395

Laser Peak Power: 25W (25,000mw)

Infrared LED's: 60mW

Red LED's: 7.5mW

Static Magnetic Field

Coverage: 4 cm²

Wavelength:

1 x Laser 905nm, 4 x Infrared 875nm, & 4 x Red 660nm

Class II Equipment Applied Part Type BF



Covers 5 times more area to treat large muscle groups

Incorporates 4 super pulsed GaAs diodes

17-9203 Extra LS50 Laser Shower Emitter

\$2,695

Laser Peak Power: 50W (50,000mw)

Infrared LED's: 60mW

Red LED's: 7.5mW

Static Magnetic Field

Coverage: 20 cm²

Wavelength:

4 x Laser 905nm, 4 x Infrared 875nm, & 4 x Red 660nm

Class II Equipment Applied Part Type BF



Highly effective and reimbursable

Combines light therapy and electrical stimulation

17-9209 Extra LaserStim Emitter

\$3,495

Laser Peak Power: 25W (25,000mw)

Infrared LED's: 60mW

Red LED's: 7.5mW

Static Magnetic Field + Electrical Stimulation

Coverage: 4 cm² (ATTENDED THERAPY ONLY)

Wavelength:

1 x Laser 905nm, 4 x Infrared 875nm, & 4 x Red 660nm

Class II Equipment Applied Part Type BF



17-9201 Set of 4 Acupuncture Point Probe

\$395

The probes allow the laser cluster to be adapted to treat smaller target areas without having to purchase additional expensive equipment.



17-9206 Basic Cart With Single Armature

\$495

Basic cart & single armature package is an ideal hands-free Therapy bundle. Our basic cart is lightweight and portable. It is just the right size for the TerraQuant and effortlessly moves from room to room on a sturdy base with wheels.

- *Lightweight, durable and easy to use*

- *Perfect for multi-tasking*



17-9207 DX Cart With Dual Arms

\$650

DX Cart with dual arms is an ultimate hands-free therapy package. DX Cart has ample storage space with a drawer and cabinet. It also has wheels, making it easy to move around the clinic.

Package includes: 2 armatures for hands-free operation , 2 C-clamps to connect armatures to the cart